

产品名称: GSTK1 Rabbit Polyclonal Antibody
产品货号: APRab11830



产品概述 (Summary)

产品名称 (Production Name)	GSTK1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GSTK1
别名 (Alternative Names)	HDCMD47P
基因 ID (Gene ID)	373156.0
蛋白 ID (SwissProt ID)	Q9Y2Q3. Synthesized peptide derived from human GSTK1 AA range: 138-188

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)

This gene encodes a member of the kappa class of the glutathione transferase superfamily of enzymes that function in

产品名称: GSTK1 Rabbit Polyclonal Antibody
产品货号: APRab11830




cellular detoxification. The encoded protein is localized to the peroxisome and catalyzes the conjugation of glutathione to a wide range of hydrophobic substrates facilitating the removal of these compounds from cells. Alternative splicing results in multiple transcript variants.[provided by RefSeq, Jan 2009],catalytic activity:RX + glutathione = HX + R-S-glutathione.,function:Significant glutathione conjugating activity is found only with the model substrate, 1-chloro-2,4-dinitrobenzene (CDNB),,similarity:Belongs to the GST superfamily. Kappa family.,subunit:Homodimer.,tissue specificity:Ubiquitous.,

研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Western blot analysis of lysates from 293T cells, primary antibody was diluted at 1:1000, 4° over night

注意事项 (Note)

For research use only .