

产品名称: GUSBP1 Rabbit Polyclonal Antibody
产品货号: APRab11861



产品概述 (Summary)

产品名称 (Production Name)	GUSBP1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GUSBP1
别名 (Alternative Names)	GUSBP1; SMA3; Putative beta-glucuronidase-like protein SMA3
基因 ID (Gene ID)	153561.0
蛋白 ID (SwissProt ID)	Q15486.The antiserum was produced against synthesized peptide derived from human GUSBL1. AA range:51-100

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	18kDa

研究背景 (Background)


产品名称: GUSBP1 Rabbit Polyclonal Antibody
产品货号: APRab11861




Binds GTP, GDP and GMP.skeletal system development,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of the lysates from HT-29 cells using GUSBL1 antibody.



Western Blot analysis of K562 cells using GUSBP1 Polyclonal Antibody

注意事项 (Note)

For research use only .