

产品名称: Gy 5 Rabbit Polyclonal Antibody
产品货号: APRab11869



产品概述 (Summary)

产品名称 (Production Name)	Gy 5 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GNG5
别名 (Alternative Names)	GNG5; GNGT5; Guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-5
基因 ID (Gene ID)	2787.0
蛋白 ID (SwissProt ID)	P63218.The antiserum was produced against synthesized peptide derived from human GNG5. AA range:10-59

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC/IF 1:200-1:1000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

产品名称: Gy 5 Rabbit Polyclonal Antibody
产品货号: APRab11869




研究背景 (Background)

G protein subunit gamma 5(GNG5) Homo sapiens G proteins are trimeric (alpha-beta-gamma) membrane-associated proteins that regulate flow of information from cell surface receptors to a variety of internal metabolic effectors. Interaction of a G protein with its activated receptor promotes exchange of GTP for GDP that is bound to the alpha subunit. The alpha-GTP complex dissociates from the beta-gamma heterodimer so that the subunits, in turn, may interact with and regulate effector molecules (Gilman, 1987 [PubMed 3113327]; summary by Ahmad et al., 1995) [PubMed 7606925]. [supplied by OMIM, Nov 2010], function: Guanine nucleotide-binding proteins (G proteins) are involved as a modulator or transducer in various transmembrane signaling systems. The beta and gamma chains are required for the GTPase activity, for replacement of GDP by GTP, and for G protein-effector interaction., similarity: Belongs to the G protein gamma family., subunit: G proteins are composed of 3 units, alpha, beta and gamma.,

研究领域 (Research Area)

Chemokine;

图片 (Image Data)



Immunofluorescence analysis of A549 cells, using GNG5 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .