

产品名称: HECD1 Rabbit Polyclonal Antibody  
产品货号: APRab11964



## 产品概述 (Summary)

|                        |                                  |
|------------------------|----------------------------------|
| 产品名称 (Production Name) | HECD1 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody       |
| 宿主 (Host)              | Rabbit                           |
| 应用 (Application)       | IHC, ICC/IF                      |
| 种属反应性 (Reactivity)     | Human, Mouse                     |

## 产品性能 (Performance)

|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                             |
| 修饰 (Modification)   | Unmodified                                                                               |
| 同种型 (Isotype)       | IgG                                                                                      |
| 克隆 (Clonality)      | Polyclonal                                                                               |
| 形式 (Form)           | Liquid                                                                                   |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.                |
| 纯化方式 (Purification) | Affinity purification                                                                    |

## 免疫原信息 (Immunogen)

|                        |                                                                      |
|------------------------|----------------------------------------------------------------------|
| 基因名 (Gene Name)        | HECTD1                                                               |
| 别名 (Alternative Names) | KIAA1131                                                             |
| 基因 ID (Gene ID)        | 25831.0                                                              |
| 蛋白 ID (SwissProt ID)   | Q9ULT8.Synthesized peptide derived from part region of human protein |

## 产品应用 (Application)

|                          |                                   |
|--------------------------|-----------------------------------|
| 稀释比 (Dilution Ratio)     | IHC 1:50-1:300, ICC/IF 1:50-1:200 |
| 蛋白分子量 (Molecular Weight) | 287kDa                            |

## 研究背景 (Background)


function: Probable E3 ubiquitin-protein ligase which accepts ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates., pathway: Protein modification; protein

**产品名称: HECD1 Rabbit Polyclonal Antibody**  
**产品货号: APRab11964**

ubiquitination.,similarity:Belongs to the UPL family. K-HECT subfamily.,similarity:Contains 1 HECT (E6AP-type E3 ubiquitin-protein ligase) domain.,similarity:Contains 1 MIB/HERC2 domain.,similarity:Contains 4 ANK repeats.,function:Probable E3 ubiquitin-protein ligase which accepts ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates.,pathway:Protein modification; protein ubiquitination.,similarity:Belongs to the UPL family. K-HECT subfamily.,similarity:Contains 1 HECT (E6AP-type E3 ubiquitin-protein ligase) domain.,similarity:Contains 1 MIB/HERC2 domain.,similarity:Contains 4 ANK repeats.,

## 研究领域 (Research Area)

## 图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Tris-EDTA,pH9.0 was used for antigen retrieval. 2 Antibody was diluted at 1:200 (4° overnight.3,Secondary antibody was diluted at 1:200 (room temperature, 45min) .

## 注意事项 (Note)

For research use only .