

产品名称: HEN2 Rabbit Polyclonal Antibody
产品货号: APRab11980



产品概述 (Summary)

产品名称 (Production Name)	HEN2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NHLH2 BHLHA34 HEN2 KIAA0490
别名 (Alternative Names)	
基因 ID (Gene ID)	4808.0
蛋白 ID (SwissProt ID)	Q02577.Synthesized peptide derived from human protein . at AA range: 30-110

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	14kDa

研究背景 (Background)

function:May serve as DNA-binding protein and may be involved in the control of cell-type determination, possibly within


产品名称: HEN2 Rabbit Polyclonal Antibody
产品货号: APRab11980



the developing nervous system.,similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,subunit:Efficient DNA binding requires dimerization with another bHLH protein.,function:May serve as DNA-binding protein and may be involved in the control of cell-type determination, possibly within the developing nervous system.,similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,subunit:Efficient DNA binding requires dimerization with another bHLH protein.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night

注意事项 (Note)

For research use only .