

产品名称: HLA-DO β Rabbit Polyclonal Antibody
产品货号: APRab12082



产品概述 (Summary)

产品名称 (Production Name)	HLA-DO β Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HLA-DOB
别名 (Alternative Names)	HLA-DOB; HLA class II histocompatibility antigen; DO beta chain; MHC class II antigen DOB
基因 ID (Gene ID)	3112.0
蛋白 ID (SwissProt ID)	P13765.The antiserum was produced against synthesized peptide derived from human HLA-DOB. AA range:1-50

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	30kDa

产品名称: HLA-DO β Rabbit Polyclonal Antibody
产品货号: AP Rab12082




研究背景 (Background)

HLA-DOB belongs to the HLA class II beta chain paralogues. This class II molecule is a heterodimer consisting of an alpha (DOA) and a beta chain (DOB), both anchored in the membrane. It is located in intracellular vesicles. DO suppresses peptide loading of MHC class II molecules by inhibiting HLA-DM. Class II molecules are expressed in antigen presenting cells (APC: B lymphocytes, dendritic cells, macrophages). The beta chain is approximately 26-28 kDa and its gene contains 6 exons. Exon one encodes the leader peptide, exons 2 and 3 encode the two extracellular domains, exon 4 encodes the transmembrane domain and exon 5 encodes the cytoplasmic tail. [provided by RefSeq, Jul 2008],function:Important modulator in the HLA class II restricted antigen presentation pathway by interaction with the HLA-DM molecule.,polymorphism:The following alleles of DOB are known: DOB*0101, DOB*0102, DOB*0103 and DOB*0104. The sequence shown is that of DOB*0101.,similarity:Belongs to the MHC class II family.,similarity:Contains 1 Ig-like C1-type (immunoglobulin-like) domain.,subunit:Heterodimer of an alpha chain (DOA) and a beta chain (DOB),.


研究领域 (Research Area)

Cell adhesion molecules (CAMs);Antigen processing and presentation;Intestinal immune network for IgA production;Type I diabetes mellitus;Asthma;Autoimmune thyroid disease;Systemic lupus erythematosus;Allograft rejection;Graft-versus-host disease;Viral myocarditis;

图片 (Image Data)




Western blot analysis of HLA-DOB Antibody. The lane on the right is blocked with the HLA-DOB peptide.



产品名称: HLA-DO β Rabbit Polyclonal Antibody
产品货号: APRab12082



Western blot analysis of the lysates from HeLa cells using HLA-DO β antibody.



Western Blot analysis of various cells using HLA-DO β Polyclonal Antibody diluted at 1: 1000

注意事项 (Note)

For research use only .