

产品名称: HLA-H Rabbit Polyclonal Antibody  
产品货号: APRab12091



## 产品概述 (Summary)

|                        |                                  |
|------------------------|----------------------------------|
| 产品名称 (Production Name) | HLA-H Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody       |
| 宿主 (Host)              | Rabbit                           |
| 应用 (Application)       | WB,ELISA                         |
| 种属反应性 (Reactivity)     | Human                            |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                      |
|------------------------|----------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | HLA-H                                                                                                                |
| 别名 (Alternative Names) | HLA-H; HLAH; Putative HLA class I histocompatibility antigen; alpha chain H; HLA-12.4; HLA-AR; MHC class I antigen H |
| 基因 ID (Gene ID)        | 3136.0                                                                                                               |
| 蛋白 ID (SwissProt ID)   | P01893.The antiserum was produced against synthesized peptide derived from human HLAH. AA range:71-120               |

## 产品应用 (Application)

|                          |                                      |
|--------------------------|--------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,ELISA 1:5000-1:20000 |
| 蛋白分子量 (Molecular Weight) | 40kDa                                |

产品名称: HLA-H Rabbit Polyclonal Antibody  
产品货号: APRab12091




## 研究背景 (Background)


This major histocompatibility complex gene represents a transcribed pseudogene, possibly derived from HLA-A. This gene displays extensive variation. [provided by RefSeq, May 2010],caution:Could be the product of a pseudogene.,function:Involved in the presentation of foreign antigens to the immune system.,similarity:Belongs to the MHC class I family.,similarity:Contains 1 Ig-like C1-type (immunoglobulin-like) domain.,subunit:Heterodimer of an alpha chain and a beta chain (beta-2-microglobulin),

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of lysates from LOVO cells, using HLAH Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HT-29 cells using HLAH antibody.

## 注意事项 (Note)

For research use only .