

产品名称: hnRNP CL1 Rabbit Polyclonal Antibody
产品货号: APRab12140



产品概述 (Summary)

产品名称 (Production Name)	hnRNP CL1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HNRNPCL1
别名 (Alternative Names)	HNRNPCL1; HNRPCL1; Heterogeneous nuclear ribonucleoprotein C-like 1; hnRNP C-like-1; hnRNP core protein C-like 1
基因 ID (Gene ID)	343069.0
蛋白 ID (SwissProt ID)	O60812.The antiserum was produced against synthesized peptide derived from human HNRCL. AA range:31-80

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	32kDa

产品名称: hnRNP CL1 Rabbit Polyclonal Antibody
产品货号: APRab12140



研究背景 (Background)

function: May play a role in nucleosome assembly by neutralizing basic proteins such as A and B core

hnRNPs., similarity: Belongs to the RRM HNRPC family. RALY subfamily., similarity: Contains 1 RRM (RNA recognition motif)


domain., subcellular location: Component of ribonucleosomes., function: May play a role in nucleosome assembly by

neutralizing basic proteins such as A and B core hnRNPs., similarity: Belongs to the RRM HNRPC family. RALY


subfamily., similarity: Contains 1 RRM (RNA recognition motif) domain., subcellular location: Component of ribonucleosomes.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of HNRCL Antibody. The lane on the right is blocked with the HNRCL peptide.



Western blot analysis of the lysates from RAW264.7 cells using HNRCL antibody.

注意事项 (Note)

For research use only .