

**产品名称: HoxD12 Rabbit Polyclonal Antibody**  
**产品货号: APRab12178**



## 产品概述 (Summary)

|                        |                                   |
|------------------------|-----------------------------------|
| 产品名称 (Production Name) | HoxD12 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody        |
| 宿主 (Host)              | Rabbit                            |
| 应用 (Application)       | WB,IHC,ICC/IF,ELISA               |
| 种属反应性 (Reactivity)     | Human,Mouse                       |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                           |
|------------------------|-----------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | HOXD12                                                                                                    |
| 别名 (Alternative Names) | HOXD12; HOX4H; Homeobox protein Hox-D12; Homeobox protein Hox-4H                                          |
| 基因 ID (Gene ID)        | 3238.0                                                                                                    |
| 蛋白 ID (SwissProt ID)   | P35452.The antiserum was produced against synthesized peptide derived from human HOXD12. AA range:191-240 |

## 产品应用 (Application)

|                          |                                                                         |
|--------------------------|-------------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000 |
| 蛋白分子量 (Molecular Weight) | 28kDa                                                                   |

## 研究背景 (Background)


产品名称: HoxD12 Rabbit Polyclonal Antibody  
产品货号: AP Rab12178




This gene belongs to the homeobox family of genes. The homeobox genes encode a highly conserved family of transcription factors that play an important role in morphogenesis in all multicellular organisms. Mammals possess four similar homeobox gene clusters, HOXA, HOXB, HOXC and HOXD, located on different chromosomes, consisting of 9 to 11 genes arranged in tandem. This gene is one of several homeobox HOXD genes located in a cluster on chromosome 2. Deletions that remove the entire HOXD gene cluster or the 5' end of this cluster have been associated with severe limb and genital abnormalities. The exact role of this gene has not been determined. [provided by RefSeq, Jul 2008],function:Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.,similarity:Belongs to the Abd-B homeobox family.,similarity:Contains 1 homeobox DNA-binding domain.,

## 研究领域 (Research Area)

## 图片 (Image Data)




Western blot analysis of lysates from COLO205 and HucEc cells, treated with serum 20% 15', using HOXD12 Antibody. The lane on the right is blocked with the synthesized peptide.




Western Blot analysis of various cells using HoxD12 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).

产品名称: HoxD12 Rabbit Polyclonal Antibody  
产品货号: APRab12178



Immunohistochemical analysis of paraffin-embedded Human breast cancer. Antibody was diluted at 1:100 (4°,overnight) . High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negative contrl (right) obtaned from antibody was pre-absorbed by immunogen peptide.



Immunohistochemical analysis of paraffin-embedded Human breast cancer. Antibody was diluted at 1:100 (4°,overnight) . High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negative contrl (right) obtaned from antibody was pre-absorbed by immunogen peptide.

### 注意事项 (Note)

For research use only .