

**产品名称: HRF Rabbit Polyclonal Antibody**  
**产品货号: APRab12201**



## 产品概述 (Summary)

|                        |                                |
|------------------------|--------------------------------|
| 产品名称 (Production Name) | HRF Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody     |
| 宿主 (Host)              | Rabbit                         |
| 应用 (Application)       | WB,IHC,ICC/IF,ELISA            |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                        |
|------------------------|--------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | TPT1                                                                                                   |
| 别名 (Alternative Names) | TPT1; Translationally-controlled tumor protein; TCTP; Fortilin; Histamine-releasing factor; HRF; p23   |
| 基因 ID (Gene ID)        | 7178.0                                                                                                 |
| 蛋白 ID (SwissProt ID)   | P13693.The antiserum was produced against synthesized peptide derived from human TCTP. AA range:71-120 |

## 产品应用 (Application)

|                          |                                                                           |
|--------------------------|---------------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:10000-1:20000 |
| 蛋白分子量 (Molecular Weight) | 24kDa                                                                     |

产品名称: HRF Rabbit Polyclonal Antibody  
产品货号: APRab12201




## 研究背景 (Background)

function:Involved in calcium binding and microtubule stabilization.,similarity:Belongs to the TCTP family.,subunit:Interacts with STEAP3.,tissue specificity:Found in several healthy and tumoral cells including erythrocytes, hepatocytes, macrophages, platelets, keratinocytes, erythroleukemia cells, gliomas, melanomas, hepatoblastomas, and lymphomas. It cannot be detected in kidney and renal cell carcinoma (RCC). Expressed in placenta and prostate.,function:Involved in calcium binding and microtubule stabilization.,similarity:Belongs to the TCTP family.,subunit:Interacts with STEAP3.,tissue specificity:Found in several healthy and tumoral cells including erythrocytes, hepatocytes, macrophages, platelets, keratinocytes, erythroleukemia cells, gliomas, melanomas, hepatoblastomas, and lymphomas. It cannot be detected in kidney and renal cell carcinoma (RCC). Expressed in placenta and prostate.,


## 研究领域 (Research Area)

Immunology

## 图片 (Image Data)




Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using TCTP Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from HepG2 and COLO205/HeLa cells, using TCTP Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: HRF Rabbit Polyclonal Antibody  
产品货号: APRab12201



Western Blot analysis of various cells using HRF Polyclonal Antibody

### 注意事项 (Note)

For research use only .