

产品名称: HSFX1 Rabbit Polyclonal Antibody
产品货号: APRab12229



产品概述 (Summary)

产品名称 (Production Name)	HSFX1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HSFX1
别名 (Alternative Names)	HSFX1; LW-1; HSF2; Heat shock transcription factor; X-linked
基因 ID (Gene ID)	100130086/100506164
蛋白 ID (SwissProt ID)	Q9UBD0.The antiserum was produced against synthesized peptide derived from human LW-1. AA range:10-59

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	47kDa

研究背景 (Background)


产品名称: HSFX1 Rabbit Polyclonal Antibody
产品货号: APRab12229




similarity:Belongs to the HSF family.,tissue specificity:Testis-specific.,similarity:Belongs to the HSF family.,tissue specificity:Testis-specific.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from HeLa cells, using LW-1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using HSFX1 Polyclonal Antibody

注意事项 (Note)

For research use only .