

产品名称: HSP 75 Rabbit Polyclonal Antibody
产品货号: APRab12232



产品概述 (Summary)

产品名称 (Production Name)	HSP 75 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TRAP1 HSP75
别名 (Alternative Names)	Heat shock protein 75 kDa, mitochondrial (HSP 75; TNFR-associated protein 1; Tumor necrosis factor type 1 receptor-associated protein; TRAP-1)
基因 ID (Gene ID)	10131.0
蛋白 ID (SwissProt ID)	Q12931. The antiserum was produced against synthesized peptide derived from the Internal region of human TRAP1. AA range: 481-530

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:200, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

产品名称: HSP 75 Rabbit Polyclonal Antibody
产品货号: APRab12232




研究背景 (Background)

This gene encodes a mitochondrial chaperone protein that is member of the heat shock protein 90 (HSP90) family. The encoded protein has ATPase activity and interacts with tumor necrosis factor type I. This protein may function in regulating cellular stress responses. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Jan 2013],function:Chaperone that expresses an ATPase activity.,similarity:Belongs to the heat shock protein 90 family.,subunit:Binds to the intracellular domain of tumor necrosis factor type 1 receptor. Binds to RB1.,tissue specificity:Found in skeletal muscle, liver, heart, brain, kidney, pancreas, lung and placenta.,


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)




Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200

产品名称: HSP 75 Rabbit Polyclonal Antibody
产品货号: APRab12232



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200

注意事项 (Note)

For research use only .