

产品名称: HSP90A Rabbit Polyclonal Antibody
产品货号: APRab12262



产品概述 (Summary)

产品名称 (Production Name)	HSP90A Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HSP90AA1
别名 (Alternative Names)	HSP90AA1; HSP90A; HSPC1; HSPCA; Heat shock protein HSP 90-alpha; Heat shock 86 kDa; HSP 86; HSP86; Renal carcinoma antigen NY-REN-38
基因 ID (Gene ID)	3320.0
蛋白 ID (SwissProt ID)	P07900.The antiserum was produced against synthesized peptide derived from human HSP90A. AA range:683-732

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	85kDa

产品名称: HSP90A Rabbit Polyclonal Antibody
产品货号: APRab12262




研究背景 (Background)

The protein encoded by this gene is an inducible molecular chaperone that functions as a homodimer. The encoded protein aids in the proper folding of specific target proteins by use of an ATPase activity that is modulated by co-chaperones. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2012],caution:Despite classification as a pseudogene, the existence of this protein is supported by unambiguous mass spectrometry evidence.,function:Molecular chaperone.,function:Molecular chaperone. Has ATPase activity.,similarity:Belongs to the heat shock protein 90 family.,subcellular location:Identified by mass spectrometry in melanosome fractions from stage I to stage IV.,subunit:Homodimer. Interacts with AHSA1, SMYD3 and TOM34. Interacts with FNIP1 and HSF1.,


研究领域 (Research Area)

PI3K/Akt; Protein_Acetylation

图片 (Image Data)




Immunofluorescence analysis of NIH/3T3 cells, using HSP90A Antibody. The picture on the right is blocked with the synthesized peptide.




Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using HSP90A Antibody. The picture on the right is blocked with the synthesized peptide.


产品名称: HSP90A Rabbit Polyclonal Antibody
产品货号: APRab12262



Western blot analysis of lysates from NIH/3T3 cells, using HSP90A Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using HSP90A Polyclonal Antibody diluted at 1: 2000



Western Blot analysis of HeLa cells using HSP90A Polyclonal Antibody diluted at 1: 2000

注意事项 (Note)

For research use only .