

产品名称: HUNK Rabbit Polyclonal Antibody
产品货号: APRab12283



产品概述 (Summary)

产品名称 (Production Name)	HUNK Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HUNK
别名 (Alternative Names)	HUNK; MAKV; Hormonally up-regulated neu tumor-associated kinase; B19; Serine/threonine-protein kinase MAK-V
基因 ID (Gene ID)	30811.0
蛋白 ID (SwissProt ID)	P57058.The antiserum was produced against synthesized peptide derived from human HUNK. AA range:251-300

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	75kDa

产品名称: HUNK Rabbit Polyclonal Antibody
产品货号: APRab12283




研究背景 (Background)

catalytic activity:ATP + a protein = ADP + a phosphoprotein.,similarity:Belongs to the protein kinase superfamily. CAMK Ser/Thr protein kinase family. SNF1 subfamily.,similarity:Contains 1 protein kinase domain.,catalytic activity:ATP + a protein = ADP + a phosphoprotein.,similarity:Belongs to the protein kinase superfamily. CAMK Ser/Thr protein kinase family. SNF1 subfamily.,similarity:Contains 1 protein kinase domain.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of the lysates from 293 cells using HUNK antibody.

注意事项 (Note)

For research use only .