

产品名称: IFIT-5 Rabbit Polyclonal Antibody
产品货号: APRab12380



产品概述 (Summary)

产品名称 (Production Name)	IFIT-5 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IFIT5 IFIT5; ISG58; RI58; Interferon-induced protein with tetratricopeptide repeats 5;
别名 (Alternative Names)	IFIT-5; Interferon-induced 58 kDa protein; Retinoic acid- and interferon-inducible 58 kDa protein; P58
基因 ID (Gene ID)	24138.0
蛋白 ID (SwissProt ID)	Q13325.The antiserum was produced against synthesized peptide derived from human IFIT5. AA range:281-330

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	50kDa

产品名称: IFIT-5 Rabbit Polyclonal Antibody
产品货号: APRab12380




研究背景 (Background)

induction:By interferons.,similarity:Belongs to the IFIT family.,similarity:Contains 8 TPR repeats.,induction:By interferons.,similarity:Belongs to the IFIT family.,similarity:Contains 8 TPR repeats.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of the lysates from HepG2 cells using IFIT5 antibody.

注意事项 (Note)

For research use only .