

产品名称: IFITM2 Rabbit Polyclonal Antibody
产品货号: APRab12382



产品概述 (Summary)

产品名称 (Production Name)	IFITM2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IFITM2
别名 (Alternative Names)	IFITM2; Interferon-induced transmembrane protein 2; Dispanin subfamily A member 2c; DSPA2c; Interferon-inducible protein 1-8D
基因 ID (Gene ID)	10581.0
蛋白 ID (SwissProt ID)	Q01629.The antiserum was produced against synthesized peptide derived from the Internal region of human IFITM2. AA range:1-50

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	15kDa

产品名称: IFITM2 Rabbit Polyclonal Antibody
产品货号: APRab12382




研究背景 (Background)

induction:By alpha and gamma interferons.,similarity:Belongs to the CD225 family.,induction:By alpha and gamma interferons.,similarity:Belongs to the CD225 family.,

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of HeLa cells using IFITM2 Polyclonal Antibody. Antibody was diluted at 1:1000. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .