

产品名称: IFITM3 Rabbit Polyclonal Antibody
产品货号: APRab12383



产品概述 (Summary)

产品名称 (Production Name)	IFITM3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IFITM3
别名 (Alternative Names)	IFITM3; Interferon-induced transmembrane protein 3; Dispanin subfamily A member 2b; DSPA2b; Interferon-inducible protein 1-8U
基因 ID (Gene ID)	10410.0
蛋白 ID (SwissProt ID)	Q01628. The antiserum was produced against synthesized peptide derived from human IFM3. AA range: 1-50

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	15kDa

产品名称: IFITM3 Rabbit Polyclonal Antibody
产品货号: APRab12383


研究背景 (Background)

The protein encoded by this gene is an interferon-induced membrane protein that helps confer immunity to influenza A H1N1 virus, West Nile virus, and dengue virus. Two transcript variants, only one of them protein-coding, have been found for this gene. Another variant encoding an N-terminally truncated isoform has been reported, but the full-length nature of this variant has not been determined. [provided by RefSeq, May 2012], induction: By alpha and gamma interferons., similarity: Belongs to the CD225 family.,


研究领域 (Research Area)

Stem Cells

图片 (Image Data)




Immunofluorescence analysis of MCF7 cells, using IFM3 Antibody. The picture on the right is blocked with the synthesized



Western blot analysis of lysates from HepG2 and HeLa cells, using IFM3 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: IFITM3 Rabbit Polyclonal Antibody
产品货号: APRab12383



Western Blot analysis of various cells using IFITM3 Polyclonal Antibody diluted at 1: 2000



Western Blot analysis of HeLa cells using IFITM3 Polyclonal Antibody diluted at 1: 2000

注意事项 (Note)

For research use only .