

产品名称: IFN14 Rabbit Polyclonal Antibody
产品货号: APRab12386



产品概述 (Summary)

产品名称 (Production Name)	IFN14 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IFNA14
别名 (Alternative Names)	
基因 ID (Gene ID)	3448.0
蛋白 ID (SwissProt ID)	P01570.Synthesized peptide derived from human protein . at AA range: 70-150

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	20kDa

研究背景 (Background)

function:Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.,similarity:Belongs to the alpha/beta interferon

产品名称: IFN14 Rabbit Polyclonal Antibody
产品货号: AP Rab12386




family, function: Produced by macrophages, IFN- α have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase, similarity: Belongs to the alpha/beta interferon family,

研究领域 (Research Area)

Cytokine-cytokine receptor interaction; Regulation of autophagy; Antigen processing and presentation; Toll_Like; RIG-I-like receptor; Cytosolic DNA-sensing pathway; Jak_STAT; Natural killer cell mediated cytotoxicity; Autoimmune thyroid disease;

图片 (Image Data)



Western blot analysis of lysates from PC12 cells, primary antibody was diluted at 1:1000, 4° over night

注意事项 (Note)

For research use only .