

产品名称: IFN- ω Rabbit Polyclonal Antibody
产品货号: APRab12406



产品概述 (Summary)

产品名称 (Production Name)	IFN- ω Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IFNW1
别名 (Alternative Names)	IFNW1; Interferon omega-1; Interferon alpha-II-1
基因 ID (Gene ID)	3467.0
蛋白 ID (SwissProt ID)	P05000.The antiserum was produced against synthesized peptide derived from the C-terminal region of human IFNW1. AA range:146-195

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	22kDa

研究背景 (Background)

产品名称: IFN- ω Rabbit Polyclonal Antibody
产品货号: APRab12406




The protein encoded by this gene is an interferon and possesses antiviral activity. The encoded protein binds to the interferon alpha/beta receptor but not to the interferon gamma receptor. This intronless gene has several pseudogenes spread throughout the genome. [provided by RefSeq, Nov 2015],similarity:Belongs to the alpha/beta interferon family.,


研究领域 (Research Area)

Cytokine-cytokine receptor interaction;RIG-I-like receptor;Jak_STAT;


图片 (Image Data)



Western blot analysis of lysate from HeLa cells, using IFNW1 Antibody.




Western Blot analysis of HeLa cells using IFN- ω Polyclonal Antibody.. Secondary antibody was diluted at 1:20000



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:100

产品名称: IFN- ω Rabbit Polyclonal Antibody
产品货号: APRab12406



Immunohistochemical analysis of paraffin-embedded human-spleen, antibody was diluted at 1:100

注意事项 (Note)

For research use only .