

产品名称: IgD Rabbit Polyclonal Antibody
产品货号: APRab12419



产品概述 (Summary)

产品名称 (Production Name)	IgD Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IGHD
别名 (Alternative Names)	Ig delta chain C region
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P01880.Synthetic peptide from human protein at AA range: 111-160

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:200, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)

similarity:Contains 3 Ig-like (immunoglobulin-like) domains.,similarity:Contains 3 Ig-like (immunoglobulin-like) domains.,


产品名称: IgD Rabbit Polyclonal Antibody
产品货号: APRab12419



研究领域 (Research Area)

Immunology

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human-tonsil, antibody was diluted at 1:200

注意事项 (Note)

For research use only .