

产品名称: IGF-I Rabbit Polyclonal Antibody
产品货号: APRab12434



产品概述 (Summary)

产品名称 (Production Name)	IGF-I Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IGF1
别名 (Alternative Names)	IBP1; IGF-IA; IGF1A; IGF1; insulin-like growth factor 1; insulin-like growth factor IA; mechano growth factor; MGF; somatomedin C
基因 ID (Gene ID)	3479.0
蛋白 ID (SwissProt ID)	P05019. The antiserum was produced against synthesized peptide derived from human IGF-I. AA range: 100-149

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	

产品名称: IGF-I Rabbit Polyclonal Antibody
产品货号: AP Rab12434




研究背景 (Background)

The protein encoded by this gene is similar to insulin in function and structure and is a member of a family of proteins involved in mediating growth and development. The encoded protein is processed from a precursor, bound by a specific receptor, and secreted. Defects in this gene are a cause of insulin-like growth factor I deficiency. Alternative splicing results in multiple transcript variants encoding different isoforms that may undergo similar processing to generate mature protein. [provided by RefSeq, Sep 2015],disease:Defects in IGF1 are the cause of insulin-like growth factor I deficiency (IGF1 deficiency) [MIM:608747]. IGF1 deficiency is an autosomal recessive disorder characterized by growth retardation, sensorineural deafness and mental retardation.,function:The insulin-like growth factors, isolated from plasma, are structurally and functionally related to insulin but have a much higher growth-promoting activity.,online information:Insulin-like growth factor 1 entry,online information:The Singapore human mutation and polymorphism database,similarity:Belongs to the insulin family.,


研究领域 (Research Area)

Oocyte meiosis;p53;mTOR;Focal adhesion;Long-term depression;Progesterone-mediated oocyte maturation;Aldosterone-regulated sodium reabsorption;Pathways in cancer;Glioma;Prostate cancer;Melanoma;Hypertrophic cardiomyopathy (HCM);Dilated cardiomyopathy;

图片 (Image Data)



Immunohistochemistry analysis of IGF-I antibody in paraffin-embedded human pancreas tissue.



产品名称: IGF-I Rabbit Polyclonal Antibody
产品货号: APRab12434



Immunohistochemical analysis of paraffin-embedded Human pancreas. Antibody was diluted at 1:100 (4°,overnight) .
High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negative contrl (right) obtaned from
antibody was pre-absorbed by immunogen peptide.

注意事项 (Note)

For research use only .