

产品名称: IGFL2 Rabbit Polyclonal Antibody
产品货号: APRab12439



产品概述 (Summary)

产品名称 (Production Name)	IGFL2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IGFL2 UNQ645/PRO1275
别名 (Alternative Names)	
基因 ID (Gene ID)	147920.0
蛋白 ID (SwissProt ID)	Q6UWQ7.Synthesized peptide derived from part region of human protein AA range: 1-50

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	13kDa

研究背景 (Background)

IGF like family member 2(IGFL2) Homo sapiens IGFL2 belongs to the insulin-like growth factor (IGF; see MIM 147440)


产品名称: IGFL2 Rabbit Polyclonal Antibody
产品货号: APRab12439



family of signaling molecules that play critical roles in cellular energy metabolism and in growth and development, especially prenatal growth (Emtage et al., 2006 [PubMed 16890402]).[supplied by OMIM, Mar 2008],similarity:Belongs to the IGFL family.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4°over night

注意事项 (Note)

For research use only .