

产品名称: IL-17B Rabbit Polyclonal Antibody
产品货号: APRab12505



产品概述 (Summary)

产品名称 (Production Name)	IL-17B Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IL17B IL20
别名 (Alternative Names)	Interleukin-17B (IL-17B; Cytokine Zcyto7; Interleukin-20; IL-20; Neuronal interleukin-17-related factor)
基因 ID (Gene ID)	27190.0
蛋白 ID (SwissProt ID)	Q9UHF5. Synthetic peptide from human protein at AA range: 1-50

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:200, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)


产品名称: IL-17B Rabbit Polyclonal Antibody
产品货号: APRab12505

The protein encoded by this gene is a T cell-derived cytokine that shares sequence similarity with IL17. This cytokine was reported to stimulate the release of TNF alpha (TNF) and IL1 beta (IL1B) from a monocytic cell line. Immunohistochemical analysis of several nerve tissues indicated that this cytokine is primarily localized to neuronal cell bodies. Alternative splicing results in multiple splice variants. [provided by RefSeq, Dec 2015],function:Stimulates the release of tumor necrosis factor alpha and IL-1-beta from the monocytic cell line THP-1.,online information:Interleukin-17 entry,similarity:Belongs to the IL-17 family.,tissue specificity:Expressed in adult pancreas, small intestine, stomach, spinal cord and testis. Less pronounced expression in prostate, colon mucosal lining, and ovary.,

研究领域 (Research Area)

Cytokine-cytokine receptor interaction;

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human-small-intestine, antibody was diluted at 1:200

注意事项 (Note)

For research use only .