

产品名称: IL-17C Rabbit Polyclonal Antibody
产品货号: APRab12506



产品概述 (Summary)

产品名称 (Production Name)	IL-17C Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IL17C UNQ561/PRO1122
别名 (Alternative Names)	Interleukin-17C (IL-17C; Cytokine CX2)
基因 ID (Gene ID)	27189.0
蛋白 ID (SwissProt ID)	Q9P0M4. Synthetic peptide from human protein at AA range: 81-130

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:200, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)

The protein encoded by this gene is a T cell-derived cytokine that shares the sequence similarity with IL17. This cytokine


产品名称: IL-17C Rabbit Polyclonal Antibody
产品货号: AP Rab12506




was reported to stimulate the release of tumor necrosis factor alpha and interleukin 1 beta from a monocytic cell line. The expression of this cytokine was found to be restricted to activated T cells. [provided by RefSeq, Jul 2008],function:Stimulates the release of tumor necrosis factor alpha and IL-1-beta from the monocytic cell line THP-1.,online information:Interleukin-17 entry,similarity:Belongs to the IL-17 family.,

研究领域 (Research Area)


图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human-colon, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-colon, antibody was diluted at 1:200



产品名称: IL-17C Rabbit Polyclonal Antibody
产品货号: APRab12506



Immunohistochemical analysis of paraffin-embedded human-pancreas, antibody was diluted at 1:200

注意事项 (Note)

For research use only .