

产品名称: IL-17D Rabbit Polyclonal Antibody
产品货号: APRab12507



产品概述 (Summary)

产品名称 (Production Name)	IL-17D Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IL17D IL27 UNQ3096/PRO21175
别名 (Alternative Names)	Interleukin-17D (IL-17D; Interleukin-27; IL-27)
基因 ID (Gene ID)	53342.0
蛋白 ID (SwissProt ID)	Q8TAD2. Synthetic peptide from human protein at AA range: 141-190

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:200, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)

The protein encoded by this gene is a cytokine that shares the sequence similarity with IL17. The treatment of endothelial


产品名称: IL-17D Rabbit Polyclonal Antibody
产品货号: AP Rab12507




cells with this cytokine has been shown to stimulate the production of other cytokines including IL6, IL8 and CSF2/ GM-CSF. The increased expression of IL8 induced by this cytokine was found to be NF-kappa B-dependent. [provided by RefSeq, Jul 2008],function:Induces expression of IL-6, IL-8, and GM-CSF from endothelial cells.,online information:Interleukin-17 entry,similarity:Belongs to the IL-17 family.,tissue specificity:Expressed preferentially in adipose, skeletal muscle and CNS.,

研究领域 (Research Area)


图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded Human-tonsil, antibody was diluted at 1:100



Immunohistochemical analysis of paraffin-embedded Human-tonsil, antibody was diluted at 1:100



产品名称: IL-17D Rabbit Polyclonal Antibody
产品货号: APRab12507



Immunohistochemical analysis of paraffin-embedded Human-brain, antibody was diluted at 1:100

注意事项 (Note)

For research use only .