

产品名称: IL-1F7 Rabbit Polyclonal Antibody
产品货号: APRab12521



产品概述 (Summary)

产品名称 (Production Name)	IL-1F7 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IL37 Interleukin-37 (FIL1 zeta; IL-1X; Interleukin-1 family member 7; IL-1F7; Interleukin-1 homolog 4; IL-1H; IL-1H4; Interleukin-1 zeta; IL-1 zeta; Interleukin-1-related protein; IL-1RP1; Interleukin-23; IL-37)
别名 (Alternative Names)	
基因 ID (Gene ID)	27178.0
蛋白 ID (SwissProt ID)	Q9NZH6. Synthetic peptide from human protein at AA range: 101-150

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:200, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

产品名称: IL-1F7 Rabbit Polyclonal Antibody
产品货号: APRab12521




研究背景 (Background)


The protein encoded by this gene is a member of the interleukin 1 cytokine family. This cytokine can bind to, and may be a ligand for interleukin 18 receptor (IL18R1/IL-1Rrp). This cytokine also binds to interleukin 18 binding protein (IL18BP), an inhibitory binding protein of interleukin 18 (IL18), and subsequently forms a complex with IL18 receptor beta subunit, and through which it inhibits the activity of IL18. This gene along with eight other interleukin 1 family genes form a cytokine gene cluster on chromosome 2. Five alternatively spliced transcript variants encoding distinct isoforms have been reported. [provided by RefSeq, Jul 2008],function: Binds to interleukin-18 receptor (IL-18R) receptor but not to IL-1 receptor. Could be a new player in the inflammatory and immune responses mediated by the IL-18/IL-18R axis.,induction: By phorbol ester (PMA) in different cell lines.,miscellaneous: The sequence shown here is derived from an EMBL/GenBank/DBJ third party annotation (TPA) entry.,online information: Interleukin-1 entry,similarity: Belongs to the IL-1 family.,tissue specificity: Isoform A, isoform B and isoform C are expressed in testis, colon, placenta, lung and lymph node. Isoform D and isoform E were found only in testis and bone marrow. Whereas only isoform A is found in brain, only isoform B in kidney and only isoform C in heart.,

研究领域 (Research Area)

图片 (Image Data)




Immunohistochemical analysis of paraffin-embedded human-small intestine, antibody was diluted at 1:200




产品名称: IL-1F7 Rabbit Polyclonal Antibody
产品货号: APRab12521



Immunohistochemical analysis of paraffin-embedded human-small-intestine, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-colon-cancer, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-colon-cancer, antibody was diluted at 1:200

注意事项 (Note)

For research use only .