

产品名称: IL-20 Rabbit Polyclonal Antibody
产品货号: APRab12532



产品概述 (Summary)

产品名称 (Production Name)	IL-20 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IL20
别名 (Alternative Names)	IL20; ZCYTO10; Interleukin-20; IL-20; Cytokine Zcyto10
基因 ID (Gene ID)	50604.0
蛋白 ID (SwissProt ID)	Q9NYY1.The antiserum was produced against synthesized peptide derived from human IL-20. AA range:12-61

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	28kDa

研究背景 (Background)

产品名称: IL-20 Rabbit Polyclonal Antibody
产品货号: APRab12532




The protein encoded by this gene is a cytokine structurally related to interleukin 10 (IL10). This cytokine has been shown to transduce its signal through signal transducer and activator of transcription 3 (STAT3) in keratinocytes. A specific receptor for this cytokine is found to be expressed in skin and upregulated dramatically in psoriatic skin, suggesting a role for this protein in epidermal function and psoriasis. [provided by RefSeq, Jul 2008],function:Cytokine that may be involved in epidermal function and psoriasis. Acts through STAT3.,similarity:Belongs to the IL-10 family.,tissue specificity:Expressed at very low levels in skin, trachea, and other tissues.,


研究领域 (Research Area)

Cytokine-cytokine receptor interaction;Jak_STAT;

图片 (Image Data)



Western blot analysis of lysate from A549 cells, using IL-20 antibody.



Western Blot analysis of A549 cells using IL-20 Polyclonal Antibody

注意事项 (Note)

For research use only .