

产品名称: IL-20 Rabbit Polyclonal Antibody
产品货号: APRab12533



产品概述 (Summary)

产品名称 (Production Name)	IL-20 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IL20
别名 (Alternative Names)	Interleukin-20 (IL-20; Cytokine Zcyto10)
基因 ID (Gene ID)	50604.0
蛋白 ID (SwissProt ID)	Q9NYY1. Synthetic peptide from human protein at AA range: 81-130

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:200, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)

The protein encoded by this gene is a cytokine structurally related to interleukin 10 (IL10). This cytokine has been shown to

产品名称: IL-20 Rabbit Polyclonal Antibody
产品货号: APRab12533




transduce its signal through signal transducer and activator of transcription 3 (STAT3) in keratinocytes. A specific receptor for this cytokine is found to be expressed in skin and upregulated dramatically in psoriatic skin, suggesting a role for this protein in epidermal function and psoriasis. [provided by RefSeq, Jul 2008],function:Cytokine that may be involved in epidermal function and psoriasis. Acts through STAT3.,similarity:Belongs to the IL-10 family.,tissue specificity:Expressed at very low levels in skin, trachea, and other tissues.,

研究领域 (Research Area)

Cytokine-cytokine receptor interaction;Jak_STAT;

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human-pancreas, antibody was diluted at 1:200

注意事项 (Note)

For research use only .