

产品名称: IL-26 Rabbit Polyclonal Antibody
产品货号: APRab12541



产品概述 (Summary)

产品名称 (Production Name)	IL-26 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IL26
别名 (Alternative Names)	IL26; AK155; Interleukin-26; IL-26; Protein AK155
基因 ID (Gene ID)	55801.0
蛋白 ID (SwissProt ID)	Q9NPH9.Synthesized peptide derived from Interleukin-26 at AA range: 21-70

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	25kDa

研究背景 (Background)

This gene was identified by its overexpression specifically in herpesvirus samimiri-transformed T cells. The encoded protein

产品名称: IL-26 Rabbit Polyclonal Antibody
产品货号: APRab12541




is a member of the IL10 family of cytokines. It is a secreted protein and may function as a homodimer. This protein is thought to contribute to the transformed phenotype of T cells after infection by herpesvirus saimiri. [provided by RefSeq, Jul 2008],induction:By Herpesvirus saimiri infection.,online information:Interleukin-26 entry,similarity:Belongs to the IL-10 family.,subunit:Homodimer.,tissue specificity:Expressed in HVS transformed T-cells but not other T-cell lines or primary stimulated T-cells.,

研究领域 (Research Area)

Cytokine-cytokine receptor interaction;Jak_STAT;

图片 (Image Data)



Western blot analysis of RAT-MUSCLE RAT-SPLEEN 293T SH-SY5Y using IL26 antibody. Antibody was diluted at 1:2000.

Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .