

产品名称: Imp3 Rabbit Polyclonal Antibody
产品货号: APRab12591



产品概述 (Summary)

产品名称 (Production Name)	Imp3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IMP3/C15orf12/MRPS4
别名 (Alternative Names)	IMP3 C15orf12 MRPS4
基因 ID (Gene ID)	55272.0
蛋白 ID (SwissProt ID)	Q9NV31.Synthetic peptide from human protein at AA range: 135-184

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	25kDa

研究背景 (Background)

This gene encodes the human homolog of the yeast Imp3 protein. The protein localizes to the nucleoli and interacts with


产品名称: Imp3 Rabbit Polyclonal Antibody
产品货号: APRab12591




the U3 snoRNP complex. The protein contains an S4 domain. [provided by RefSeq, Jul 2008],function:Component of the 60-80S U3 small nucleolar ribonucleoprotein (U3 snoRNP). Required for the early cleavages during pre-18S ribosomal RNA processing.,similarity:Belongs to the ribosomal protein S4P family.,similarity:Contains 1 S4 RNA-binding domain.,subunit:Component of a heterotrimeric complex containing IMP3, IMP4 and MPHOSPH10. Interacts with MPHOSPH10.,

研究领域 (Research Area)

图片 (Image Data)




Western Blot analysis of HEPG2 cells using Antibody diluted at 1000. Secondary antibody was diluted at 1:20000




Immunohistochemical analysis of paraffin-embedded human-lung-cancer, antibody was diluted at 1:200


产品名称: Imp3 Rabbit Polyclonal Antibody
产品货号: APRab12591



Immunohistochemical analysis of paraffin-embedded human-lung-cancer, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200

注意事项 (Note)

For research use only .