

产品名称: INSL5 Rabbit Polyclonal Antibody
产品货号: APRab12626



产品概述 (Summary)

产品名称 (Production Name)	INSL5 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	INSL5 UNQ156/PRO182
别名 (Alternative Names)	
基因 ID (Gene ID)	10022.0
蛋白 ID (SwissProt ID)	Q9Y5Q6.Synthesized peptide derived from human protein . at AA range: 30-110

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	14kDa

研究背景 (Background)

The protein encoded by this gene contains a classical signature of the insulin superfamily and is highly similar to relaxin 3


产品名称: INSL5 Rabbit Polyclonal Antibody
产品货号: APRab12626



(RLN3/INSL7). [provided by RefSeq, Jul 2008],function:May have a role in gut contractility or in thymic development and regulation.,similarity:Belongs to the insulin family.,subunit:Heterodimer of a B chain and an A chain linked by two disulfide bonds.,tissue specificity:Highly expressed in rectum with lower levels in uterus and ascending and descending colon.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night

注意事项 (Note)

For research use only .