

产品名称: IQGAP3 Rabbit Polyclonal Antibody
产品货号: APRab12723



产品概述 (Summary)

产品名称 (Production Name)	IQGAP3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IQGAP3
别名 (Alternative Names)	IQGAP3; Ras GTPase-activating-like protein IQGAP3
基因 ID (Gene ID)	128239.0
蛋白 ID (SwissProt ID)	Q86VI3.The antiserum was produced against synthesized peptide derived from human IQGAP3. AA range:107-156

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	180kDa

研究背景 (Background)

产品名称: IQGAP3 Rabbit Polyclonal Antibody
产品货号: APRab12723




similarity:Contains 1 CH (calponin-homology) domain.,similarity:Contains 1 Ras-GAP domain.,similarity:Contains 4 IQ domains.,similarity:Contains 1 CH (calponin-homology) domain.,similarity:Contains 1 Ras-GAP domain.,similarity:Contains 4 IQ domains.,


研究领域 (Research Area)

Regulates Actin and Cytoskeleton;

图片 (Image Data)



Western blot analysis of the lysates from 293 cells using IQGAP3 antibody.



Western Blot analysis of various cells using IQGAP3 Polyclonal Antibody diluted at 1: 1000

注意事项 (Note)

For research use only .