

产品名称: Jamip2 Rabbit Polyclonal Antibody  
产品货号: AP Rab12828



## 产品概述 (Summary)

|                        |                                   |
|------------------------|-----------------------------------|
| 产品名称 (Production Name) | Jamip2 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody        |
| 宿主 (Host)              | Rabbit                            |
| 应用 (Application)       | WB,IHC,ICC/IF,ELISA               |
| 种属反应性 (Reactivity)     | Human,Rat,Mouse                   |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                                                                         |
|------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | JAKMIP2<br>JAKMIP2; JAMIP2; KIAA0555; NECC1; Janus kinase and microtubule-interacting protein 2; CTCL tumor antigen HD-CL-04; Neuroendocrine long coiled-coil protein 1 |
| 别名 (Alternative Names) |                                                                                                                                                                         |
| 基因 ID (Gene ID)        | 9832.0                                                                                                                                                                  |
| 蛋白 ID (SwissProt ID)   | Q96AA8.The antiserum was produced against synthesized peptide derived from human JAKMIP2. AA range:761-810                                                              |

## 产品应用 (Application)

|                          |                                                                          |
|--------------------------|--------------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000 |
| 蛋白分子量 (Molecular Weight) | 95kDa                                                                    |

产品名称: Jamip2 Rabbit Polyclonal Antibody  
产品货号: APRab12828




## 研究背景 (Background)


The protein encoded by this gene is reported to be a component of the Golgi matrix. It may act as a golgin protein by negatively regulating transit of secretory cargo and by acting as a structural scaffold of the Golgi. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2012],similarity:Belongs to the JAKMIP family.,tissue specificity:Highly expressed in brain, moderately expressed in thymus, spleen and lung, and weakly expressed in kidney, liver and peripheral blood lymphocytes. Also expressed in adrenal and pituitary glands, as well as testis.,

## 研究领域 (Research Area)

## 图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using JAKMIP2 Antibody. The picture on the right is blocked with the synthesized peptide.



## 注意事项 (Note)

For research use only .

产品名称: Jamip2 Rabbit Polyclonal Antibody  
产品货号: APRab12828

---

