

产品名称: JTB Rabbit Polyclonal Antibody
产品货号: APRab12851



产品概述 (Summary)

产品名称 (Production Name)	JTB Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	JTB
别名 (Alternative Names)	JTB; HSPC222; Protein JTB; Jumping translocation breakpoint protein; Prostate androgen-regulated protein; PAR protein
基因 ID (Gene ID)	10899.0
蛋白 ID (SwissProt ID)	O76095.The antiserum was produced against synthesized peptide derived from human JTB. AA range:10-59

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	17kDa

产品名称: JTB Rabbit Polyclonal Antibody
产品货号: APRab12851




研究背景 (Background)


similarity:Belongs to the JTB family.,tissue specificity:Expressed in all normal human tissues studied but overexpressed in most of their malignant counterparts.,similarity:Belongs to the JTB family.,tissue specificity:Expressed in all normal human tissues studied but overexpressed in most of their malignant counterparts.,

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of JTB Antibody. The lane on the right is blocked with the JTB peptide.



Western blot analysis of the lysates from HUVEC cells using JTB antibody.



产品名称: JTB Rabbit Polyclonal Antibody
产品货号: APRab12851



Western Blot analysis of k562 cells using JTB Polyclonal Antibody

注意事项 (Note)

For research use only .