

产品名称: Jun B Rabbit Polyclonal Antibody
产品货号: APRab12852



产品概述 (Summary)

产品名称 (Production Name)	Jun B Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	JUNB
别名 (Alternative Names)	JUNB; Transcription factor jun-B
基因 ID (Gene ID)	3726.0
蛋白 ID (SwissProt ID)	P17275.The antiserum was produced against synthesized peptide derived from human JunB. AA range:226-275

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	36kDa

产品名称: Jun B Rabbit Polyclonal Antibody
产品货号: APRab12852




研究背景 (Background)

function: Transcription factor involved in regulating gene activity following the primary growth factor response. Binds to the DNA sequence 5'-TGA[CG]TCA-3', induction: By growth factors., similarity: Belongs to the bZIP family., similarity: Belongs to the bZIP family. Jun subfamily., similarity: Contains 1 bZIP domain., subunit: Binds DNA as an homodimer or as an heterodimer with another member of the jun/fos family., function: Transcription factor involved in regulating gene activity following the primary growth factor response. Binds to the DNA sequence 5'-TGA[CG]TCA-3', induction: By growth factors., similarity: Belongs to the bZIP family., similarity: Belongs to the bZIP family. Jun subfamily., similarity: Contains 1 bZIP domain., subunit: Binds DNA as an homodimer or as an heterodimer with another member of the jun/fos family.,


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)




Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using JunB Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from HeLa cells, using JunB Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: Jun B Rabbit Polyclonal Antibody
产品货号: APRab12852



Western Blot analysis of HeLa RAW264.7 cells using Jun B Polyclonal Antibody cells nucleus extracted by Minute TM
Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA) .

注意事项 (Note)

For research use only .