

**产品名称:** KI2LA Rabbit Polyclonal Antibody

**产品货号:** APRab12997

### 产品概述 (Summary)

|                   |                                                                                          |
|-------------------|------------------------------------------------------------------------------------------|
| 描述(Description)   | Rabbit polyclonal Antibody                                                               |
| 宿主(Host)          | Rabbit                                                                                   |
| 应用(Application)   | WB,ELISA                                                                                 |
| 种属反应(Reactivity)  | Human,Rat,Mouse                                                                          |
| 偶联物(Conjugation)  | Unconjugated                                                                             |
| 修饰(Modification)  | Unmodified                                                                               |
| 同种型(Isotype)      | IgG                                                                                      |
| 克隆(Clonality)     | Polyclonal                                                                               |
| 剂型(Form)          | Liquid                                                                                   |
| 保存条件(Storage)     | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| 储存溶液(Buffer)      | Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.                |
| 纯化 (Purification) | Affinity purification                                                                    |

### 产品应用(Application)

|                       |                                      |
|-----------------------|--------------------------------------|
| 稀释比(Dilution Ratio)   | WB 1:500-1:2000,ELISA 1:5000-1:20000 |
| 分子量(Molecular Weight) | 41kDa                                |

### 抗原信息(Antigen Information)

|                       |                                                               |
|-----------------------|---------------------------------------------------------------|
| 基因名(Gene Name)        | KIR2DL5A CD158F CD158F1 KIR2DL5                               |
| 别名(Alternative Names) |                                                               |
| 基因 ID(Gene ID)        | 57292.0                                                       |
| SwissProt ID          | Q8N109                                                        |
| 免疫原(Immunogen)        | Synthesized peptide derived from part region of human protein |

### 研究背景 (Background)

killer cell immunoglobulin like receptor, two Ig domains and long cytoplasmic tail 5A(KIR2DL5A) Homo sapiens Killer cell immunoglobulin-like receptors (KIRs) are transmembrane glycoproteins expressed by natural killer cells and subsets of T cells. The KIR genes are polymorphic and highly homologous and they are found in a cluster on chromosome 19q13.4 within the 1 Mb leukocyte receptor complex (LRC). The gene content of the KIR gene cluster varies among haplotypes, although several

"framework" genes are found in all haplotypes (KIR3DL3, KIR3DP1, KIR3DL4, KIR3DL2). The KIR proteins are classified by the number of extracellular immunoglobulin domains (2D or 3D) and by whether they have a long (L) or short (S) cytoplasmic domain. KIR proteins with the long cytoplasmic domain transduce inhibitory signals upon ligand binding via an immune tyrosine-based inhibitory motif (ITIM), while KIR proteins with the short cytoplasmic domain lack the function: Receptor on natural killer (NK) cells for HLA-C alleles. Inhibits the activity of NK cells thus preventing cell lysis., similarity: Belongs to the immunoglobulin superfamily., similarity: Contains 2 Ig-like C2-type (immunoglobulin-like) domains.,


### 研究领域 (Research Area)

Antigen processing and presentation; Natural killer cell mediated cytotoxicity; Graft-versus-host disease;

### 注意事项 (Note)

For research use only.

### 图片 (Image Data)



Western blot analysis of lysates from PC12 cells, KI2LA Rabbit Polyclonal Antibody was diluted at 1:1000, 4° over night

## EnkiLife 优势产品推荐

| 产品类别                        | 产品货号     | 产品名称                                               |
|-----------------------------|----------|----------------------------------------------------|
| WB 解决方案                     | RA10020  | 2h 极速 WB 即用型全流程试剂盒                                 |
|                             | RA10021  | 4h 快速 WB 即用型全流程试剂盒                                 |
|                             | RA10042  | 5h 畅享版 WB 全流程试剂盒                                   |
|                             | RA10037  | 校准级彩色预染蛋白 Marker (8-180kDa)                        |
|                             | RA10038  | 校准级彩色预染蛋白 Marker (10-250kDa)                       |
|                             | RA10039  | 校准级高分子彩色预染蛋白 Marker (25-400kDa)                    |
| TSA 多重荧光染色试剂盒               | RA10008  | TSA 双标三色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10009  | TSA 三标四色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10010  | TSA 四标五色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10011  | TSA 五标六色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10012  | TSA 六标七色多重荧光染色试剂盒 (mIHC)                           |
| IHC 检测试剂盒                   | RA10006  | HRP Anti-Mouse/Rabbit IHC Detection System         |
|                             | RA10007  | Polymer-HRP Anti-Mouse/Rabbit IHC Detection System |
| 抗体标记试剂盒                     | RE80004p | 辣根过氧化物酶(HRP)抗体标记试剂盒                                |
|                             | RE80002q | Sulfo-NHS-生物素标记试剂盒                                 |
|                             | RE80007p | Cy3 荧光素标记试剂盒                                       |
|                             | RE80011p | Fluor488 荧光素标记试剂盒                                  |
|                             | RE80017p | Fluor750 荧光素标记试剂盒                                  |
|                             | RE80005p | 藻红蛋白(R-PE) 抗体快速标记试剂盒                               |
|                             | RE80040  | PE-Cy7 串联染料抗体快速标记试剂盒                               |
| 稳转细胞系构建服务<br>(免费赠送全膜 WB 验证) | TS-0001  | 过表达稳转细胞系构建                                         |
|                             | TS-0002  | 敲低稳转细胞系构建                                          |
|                             | TS-0003  | 敲除细胞系构建                                            |

限于篇幅，更多产品详见官网：<https://www.enkilife.cn>