

产品名称: KLHL29 Rabbit Polyclonal Antibody
产品货号: APRab13068



产品概述 (Summary)

产品名称 (Production Name)	KLHL29 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KLHL29
别名 (Alternative Names)	KLHL29; KBTBD9; KIAA1921; Kelch-like protein 29; Kelch repeat and BTB domain-containing protein 9
基因 ID (Gene ID)	114818.0
蛋白 ID (SwissProt ID)	Q96CT2. The antiserum was produced against synthesized peptide derived from human KLHL29. AA range: 261-310

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	

产品名称: KLHL29 Rabbit Polyclonal Antibody
产品货号: APRab13068




研究背景 (Background)

caution:Although the complete sequence is not known with certainty, sequence shown here appears to be the most probable in accordance with the mouse sequence ortholog.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 6 Kelch repeats.,caution:Although the complete sequence is not known with certainty, sequence shown here appears to be the most probable in accordance with the mouse sequence ortholog.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 6 Kelch repeats.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using KLHL29 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .