

产品名称: KOX4 Rabbit Polyclonal Antibody
产品货号: APRab13093



产品概述 (Summary)

产品名称 (Production Name)	KOX4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZNF7
别名 (Alternative Names)	ZNF7; KOX4; Zinc finger protein 7; Zinc finger protein HF.16; Zinc finger protein KOX4
基因 ID (Gene ID)	7553.0
蛋白 ID (SwissProt ID)	P17097. The antiserum was produced against synthesized peptide derived from human ZNF7. AA range: 41-90

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:200-1:1000, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

产品名称: KOX4 Rabbit Polyclonal Antibody
产品货号: APRab13093




研究背景 (Background)


function: May be involved in transcriptional regulation., similarity: Belongs to the krueppel C2H2-type zinc-finger protein family., similarity: Contains 1 KRAB domain., similarity: Contains 15 C2H2-type zinc fingers., tissue specificity: Ubiquitously present in many human cell lines of different embryological derivation., function: May be involved in transcriptional regulation., similarity: Belongs to the krueppel C2H2-type zinc-finger protein family., similarity: Contains 1 KRAB domain., similarity: Contains 15 C2H2-type zinc fingers., tissue specificity: Ubiquitously present in many human cell lines of different embryological derivation.,

研究领域 (Research Area)

图片 (Image Data)



Immunofluorescence analysis of HeLa cells, using ZNF7 Antibody. The picture on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using ZNF7 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .

产品名称: KOX4 Rabbit Polyclonal Antibody
产品货号: APRab13093

