

产品名称: Laminin γ -1 Rabbit Polyclonal Antibody
产品货号: AP Rab13203



产品概述 (Summary)

产品名称 (Production Name)	Laminin γ -1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat,Monkey,Cat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LAMC1
别名 (Alternative Names)	LAMC1; LAMB2; Laminin subunit gamma-1; Laminin B2 chain; Laminin-1 subunit gamma; Laminin-10 subunit gamma; Laminin-11 subunit gamma; Laminin-2 subunit gamma; Laminin-3 subunit gamma; Laminin-4 subunit gamma; Laminin-6 subunit gamma; Lamini
基因 ID (Gene ID)	3915.0
蛋白 ID (SwissProt ID)	P11047.The antiserum was produced against synthesized peptide derived from human LAMC1. AA range:1451-1500

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:20000-1:40000
----------------------	---

产品名称: Laminin γ -1 Rabbit Polyclonal Antibody
产品货号: AP Rab13203



蛋白分子量 (Molecular Weight) 178kDa


研究背景 (Background)

Laminins, a family of extracellular matrix glycoproteins, are the major noncollagenous constituent of basement membranes. They have been implicated in a wide variety of biological processes including cell adhesion, differentiation, migration, signaling, neurite outgrowth and metastasis. Laminins, composed of 3 non identical chains: laminin alpha, beta and gamma (formerly A, B1, and B2, respectively), have a cruciform structure consisting of 3 short arms, each formed by a different chain, and a long arm composed of all 3 chains. Each laminin chain is a multidomain protein encoded by a distinct gene. Several isoforms of each chain have been described. Different alpha, beta and gamma chain isomers combine to give rise to different heterotrimeric laminin isoforms which are designated by Arabic numerals in the order of their discovery, i.e. $\alpha 1\beta 1\gamma 1$ heterotrimer is laminin 1. The biological function: Domains VI and IV are globular, domain: The alpha-helical domains I and II are thought to interact with other laminin chains to form a coiled coil structure, function: Binding to cells via a high affinity receptor, laminin is thought to mediate the attachment, migration and organization of cells into tissues during embryonic development by interacting with other extracellular matrix components, similarity: Contains 1 laminin IV type A domain, similarity: Contains 1 laminin N-terminal domain, similarity: Contains 11 laminin EGF-like domains, subunit: Laminin is a complex glycoprotein, consisting of three different polypeptide chains (alpha, beta, gamma), which are bound to each other by disulfide bonds into a cross-shaped molecule comprising one long and three short arms with globules at each end. Gamma-1 is a subunit of laminin-1 (EHS laminin), laminin-2 (merosin), laminin-3 (S-laminin), laminin-4 (S-merosin), laminin-6 (K-laminin) and laminin-7 (KS-laminin), tissue specificity: Found in the basement membranes (major component),

研究领域 (Research Area)


Focal adhesion; ECM-receptor interaction; Prion diseases; Pathways in cancer; Small cell lung cancer;

图片 (Image Data)



Western blot analysis of lysates from HUVEC cells, using LAMC1 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: Laminin γ -1 Rabbit Polyclonal Antibody
产品货号: APRab13203



Western Blot analysis of customer's (cat sample) using Laminin γ -1 Polyclonal Antibody diluted at 1: 1000

注意事项 (Note)

For research use only .