

产品名称: LDLRAD2 Rabbit Polyclonal Antibody
产品货号: APRab13270



产品概述 (Summary)

产品名称 (Production Name)	LDLRAD2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LDLRAD2
别名 (Alternative Names)	LDLRAD2; Low-density lipoprotein receptor class A domain-containing protein 2
基因 ID (Gene ID)	401944.0
蛋白 ID (SwissProt ID)	Q5SZI1.The antiserum was produced against synthesized peptide derived from human LDLRAD2. AA range:130-179

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	29kDa

产品名称: LDLRAD2 Rabbit Polyclonal Antibody
产品货号: APRab13270




研究背景 (Background)


similarity:Belongs to the LDLR family.,similarity:Contains 1 LDL-receptor class A domain.,similarity:Belongs to the LDLR family.,similarity:Contains 1 LDL-receptor class A domain.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from Jurkat, LOVO, and NIH/3T3 cells, using LDLRAD2 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HepG2 cells using LDLRAD2 antibody.



产品名称: LDLRAD2 Rabbit Polyclonal Antibody
产品货号: APRab13270



Western Blot analysis of various cells using LDLRAD2 Polyclonal Antibody diluted at 1: 500

注意事项 (Note)

For research use only .