

产品名称: LDOC1 Rabbit Polyclonal Antibody
产品货号: APRab13272



产品概述 (Summary)

产品名称 (Production Name)	LDOC1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LDOC1
别名 (Alternative Names)	LDOC1; BCUR1; Protein LDOC1; Leucine zipper protein down-regulated in cancer cells
基因 ID (Gene ID)	23641.0
蛋白 ID (SwissProt ID)	O95751. The antiserum was produced against synthesized peptide derived from human LDOC1. AA range: 21-70

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	

产品名称: LDOC1 Rabbit Polyclonal Antibody
产品货号: APRab13272




研究背景 (Background)

The protein encoded by this gene contains a leucine zipper-like motif and a proline-rich region that shares marked similarity with an SH3-binding domain. The protein localizes to the nucleus and is down-regulated in some cancer cell lines. It is thought to regulate the transcriptional response mediated by the nuclear factor kappa B (NF-kappaB). The gene has been proposed as a tumor suppressor gene whose protein product may have an important role in the development and/or progression of some cancers. [provided by RefSeq, Jul 2008],function:May have an important role in the development and/or progression of some cancers.,similarity:Belongs to the LDOC1 family.,tissue specificity:Ubiquitously expressed with high levels in brain ant thyroid and low expression in placenta, liver and leukocytes. Expressed as well in six of the seven human breast cancer cell lines examined.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human Colon cancer. 1, Antibody was diluted at 1:200 (4° overnight) .
2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200 (room temperature, 45min) .

注意事项 (Note)

For research use only .