

产品名称: LDOC1L Rabbit Polyclonal Antibody  
产品货号: APRab13273



## 产品概述 (Summary)

|                        |                                   |
|------------------------|-----------------------------------|
| 产品名称 (Production Name) | LDOC1L Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody        |
| 宿主 (Host)              | Rabbit                            |
| 应用 (Application)       | WB,IHC,ICC/IF,ELISA               |
| 种属反应性 (Reactivity)     | Human,Mouse                       |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                                              |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | LDOC1L                                                                                                                                       |
| 别名 (Alternative Names) | LDOC1L; MAR6; MART6; Protein LDOC1L; Leucine zipper protein down-regulated in cancer cells-like; Mammalian retrotransposon-derived protein 6 |
| 基因 ID (Gene ID)        | 84247.0                                                                                                                                      |
| 蛋白 ID (SwissProt ID)   | Q6ICC9.The antiserum was produced against synthesized peptide derived from human LDOC1L. AA range:141-190                                    |

## 产品应用 (Application)

|                          |                                                                         |
|--------------------------|-------------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000 |
| 蛋白分子量 (Molecular Weight) | 24kDa                                                                   |

产品名称: LDOC1L Rabbit Polyclonal Antibody  
产品货号: APRab13273




## 研究背景 (Background)


similarity:Belongs to the LDOC1 family.,similarity:Belongs to the LDOC1 family.,

## 研究领域 (Research Area)


## 图片 (Image Data)



Western blot analysis of lysates from COLO205 cells, using LDOC1L Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using LDOC1L Polyclonal Antibody



Immunohistochemical analysis of paraffin-embedded Human brain. Antibody was diluted at 1:100 (4°,overnight) . High-

**产品名称: LDOC1L Rabbit Polyclonal Antibody**  
**产品货号: APRab13273**



---

pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negative control (right) obtained from antibody was pre-absorbed by immunogen peptide.

### **注意事项 (Note)**

For research use only .