

产品名称: LIR-6 Rabbit Polyclonal Antibody
产品货号: AP Rab13333



产品概述 (Summary)

产品名称 (Production Name)	LIR-6 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LILRA1 LILRA1; LIR6; Leukocyte immunoglobulin-like receptor subfamily A member 1;
别名 (Alternative Names)	CD85 antigen-like family member I; Leukocyte immunoglobulin-like receptor 6; LIR-6; CD antigen CD85i
基因 ID (Gene ID)	11024.0
蛋白 ID (SwissProt ID)	O75019.The antiserum was produced against synthesized peptide derived from human LILRA1. AA range:53-102

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	53kDa

产品名称: LIR-6 Rabbit Polyclonal Antibody
产品货号: APRab13333




研究背景 (Background)


This gene encodes an activating member of the leukocyte immunoglobulin-like receptor (LIR) family, which is found in a gene cluster at chromosomal region 19q13.4. The encoded protein is predominantly expressed in B cells, interacts with major histocompatibility complex class I ligands, and contributes to the regulation of immune responses. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, May 2013],function:May act as receptor for class I MHC antigens.,similarity:Contains 4 Ig-like C2-type (immunoglobulin-like) domains.,tissue specificity:Detected in monocytes and B-cells.,

研究领域 (Research Area)

图片 (Image Data)




Immunofluorescence analysis of A549 cells, using LILRA1 Antibody. The picture on the right is blocked with the synthesized




Western blot analysis of lysates from HeLa and MCF-7 cells, using LILRA1 Antibody. The lane on the right is blocked with the synthesized peptide.


产品名称: LIR-6 Rabbit Polyclonal Antibody
产品货号: APRab13333



Western blot analysis of the lysates from HeLa cells using LILRA1 antibody.



Western Blot analysis of various cells using LIR-6 Polyclonal Antibody



Western Blot analysis of MCF7 cells using LIR-6 Polyclonal Antibody

注意事项 (Note)

For research use only .