

产品名称: LMO4 Rabbit Polyclonal Antibody
产品货号: APRab13359



产品概述 (Summary)

产品名称 (Production Name)	LMO4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LMO4
别名 (Alternative Names)	LMO4; LIM domain transcription factor LMO4; Breast tumor autoantigen; LIM domain only protein 4; LMO-4
基因 ID (Gene ID)	8543.0
蛋白 ID (SwissProt ID)	P61968.The antiserum was produced against synthesized peptide derived from human LMO4. AA range:116-165

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	44kDa

产品名称: LMO4 Rabbit Polyclonal Antibody
产品货号: APRab13359




研究背景 (Background)


LIM domain only 4(LMO4) Homo sapiens This gene encodes a cysteine-rich protein that contains two LIM domains but lacks a DNA-binding homeodomain. The encoded protein may play a role as a transcriptional regulator or as an oncogene. [provided by RefSeq, Aug 2008],function:Probable transcriptional factor.,similarity:Contains 2 LIM zinc-binding domains.,subunit:Interacts strongly with LDBS. Interacts with DEAF1, CLIM1 and CLIM2.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from 293 cells, using LMO4 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HepG2 cells using LMO4 antibody.

产品名称: LMO4 Rabbit Polyclonal Antibody
产品货号: APRab13359



Western Blot analysis of various cells using LMO4 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA) .

注意事项 (Note)

For research use only .