

产品名称: LMTK3 Rabbit Polyclonal Antibody
产品货号: APRab13363



产品概述 (Summary)

产品名称 (Production Name)	LMTK3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LMTK3
别名 (Alternative Names)	LMTK3; KIAA1883; TYKLM3; Serine/threonine-protein kinase LMTK3; Lemur tyrosine kinase 3
基因 ID (Gene ID)	114783.0
蛋白 ID (SwissProt ID)	Q96Q04.The antiserum was produced against synthesized peptide derived from human LMTK3. AA range:1251-1300

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	40kDa

产品名称: LMTK3 Rabbit Polyclonal Antibody
产品货号: APRab13363




研究背景 (Background)


catalytic activity:ATP + a protein = ADP + a phosphoprotein.,similarity:Belongs to the protein kinase superfamily. Tyr protein kinase family.,similarity:Contains 1 protein kinase domain.,catalytic activity:ATP + a protein = ADP + a phosphoprotein.,similarity:Belongs to the protein kinase superfamily. Tyr protein kinase family.,similarity:Contains 1 protein kinase domain.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from 293 cells, using LMTK3 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HeLa cells using LMTK3 antibody.

产品名称: LMTK3 Rabbit Polyclonal Antibody
产品货号: APRab13363



Western Blot analysis of 3T3 cells using LMTK3 Polyclonal Antibody diluted at 1: 2000

注意事项 (Note)

For research use only .