

产品名称: LOST1 Rabbit Polyclonal Antibody
产品货号: APRab13373



产品概述 (Summary)

产品名称 (Production Name)	LOST1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TUSC5 TUSC5; IFITMD3; LOST1; Tumor suppressor candidate 5; Dispanin subfamily B member 1; DSPB1; Interferon-induced transmembrane domain-containing protein D3; Protein located at seventeen-p-thirteen point three 1
别名 (Alternative Names)	
基因 ID (Gene ID)	286753.0
蛋白 ID (SwissProt ID)	Q8IXB3.The antiserum was produced against synthesized peptide derived from human TUSC5. AA range:1-50

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	22kDa

产品名称: LOST1 Rabbit Polyclonal Antibody
产品货号: APRab13373




研究背景 (Background)


developmental stage: Expressed in fetal brain., function: May be involved in fat metabolism., similarity: Belongs to the CD225 family., tissue specificity: Expressed at high levels in heart, mammary gland, adrenal gland, stomach, smooth muscle and skeletal muscle, and at lower levels in brain and lung. Strongly down-regulated in lung cancer tissues, due to hypermethylation of the corresponding locus., developmental stage: Expressed in fetal brain., function: May be involved in fat metabolism., similarity: Belongs to the CD225 family., tissue specificity: Expressed at high levels in heart, mammary gland, adrenal gland, stomach, smooth muscle and skeletal muscle, and at lower levels in brain and lung. Strongly down-regulated in lung cancer tissues, due to hypermethylation of the corresponding locus.,

研究领域 (Research Area)

图片 (Image Data)



Immunofluorescence analysis of HeLa cells, using TUSC5 Antibody. The picture on the right is blocked with the synthesized



Western blot analysis of lysates from LOVO cells, using TUSC5 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: LOST1 Rabbit Polyclonal Antibody
产品货号: APRab13373



注意事项 (Note)

For research use only .