

产品名称: LRG1 Rabbit Polyclonal Antibody
产品货号: APRab13419



产品概述 (Summary)

产品名称 (Production Name)	LRG1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LRG1
别名 (Alternative Names)	LRG1; LRG; Leucine-rich alpha-2-glycoprotein; LRG
基因 ID (Gene ID)	116844.0
蛋白 ID (SwissProt ID)	P02750.Synthesized peptide derived from the C-terminal region of human LRG1.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	58kDa

研究背景 (Background)

产品名称: LRG1 Rabbit Polyclonal Antibody
产品货号: AP Rab13419




The leucine-rich repeat (LRR) family of proteins, including LRG1, have been shown to be involved in protein-protein interaction, signal transduction, and cell adhesion and development. LRG1 is expressed during granulocyte differentiation (O'Donnell et al., 2002 [PubMed 12223515]).[supplied by OMIM, Mar 2008],similarity:Contains 8 LRR (leucine-rich) repeats.,tissue specificity:Plasma.,


研究领域 (Research Area)

Immunology

图片 (Image Data)



Western Blot analysis of various cells using LRG1 Polyclonal Antibody



Western Blot analysis of HeLa cells using LRG1 Polyclonal Antibody

注意事项 (Note)

For research use only .