

产品名称: M6A Rabbit Polyclonal Antibody
产品货号: APRab13538



产品概述 (Summary)

产品名称 (Production Name)	M6A Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPM6A
别名 (Alternative Names)	GPM6A; M6A; Neuronal membrane glycoprotein M6-a; M6a
基因 ID (Gene ID)	2823.0
蛋白 ID (SwissProt ID)	P51674.Synthesized peptide derived from M6A . at AA range: 140-220

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	31kDa

研究背景 (Background)

similarity:Belongs to the myelin proteolipid protein family.,similarity:Belongs to the myelin proteolipid protein family.,


产品名称: M6A Rabbit Polyclonal Antibody
产品货号: AP Rab13538



研究领域 (Research Area)

Metabolism; Types of disease; Obesity; Epigenetics and Nuclear Signaling; RNA; RNA / Nucleotides

图片 (Image Data)



Western Blot analysis of extracts from rat stomach, using M6A Polyclonal Antibody.. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .