

产品名称: MAGC3 Rabbit Polyclonal Antibody
产品货号: APRab13584



产品概述 (Summary)

产品名称 (Production Name)	MAGC3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MAGEC3 HCA2
别名 (Alternative Names)	
基因 ID (Gene ID)	139081.0
蛋白 ID (SwissProt ID)	Q8TD91.Synthesized peptide derived from human protein . at AA range: 60-140

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	70kDa

研究背景 (Background)

This gene is a member of the MAGEC gene family. The members of this family are not expressed in normal tissues, except


产品名称: MAGC3 Rabbit Polyclonal Antibody
产品货号: AP Rab13584



for testis, and are expressed in tumors of various histological types. The MAGEC genes are clustered on chromosome Xq26-q27. Two transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jul 2008], similarity: Contains 2 MAGE domains., tissue specificity: Expressed in testis. Not expressed in other normal tissues, but is expressed in tumors of different histological origins.,

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of HEK293 lysis, using primary antibody at 1:1000 dilution. Secondary antibody was diluted at 1:10000

注意事项 (Note)

For research use only .