

产品名称: MAGH1 Rabbit Polyclonal Antibody
产品货号: APRab13598



产品概述 (Summary)

产品名称 (Production Name)	MAGH1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MAGEH1 APR1
别名 (Alternative Names)	
基因 ID (Gene ID)	28986.0
蛋白 ID (SwissProt ID)	Q9H213.Synthesized peptide derived from human protein . at AA range: 10-90

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	24kDa

研究背景 (Background)

This gene belongs to the non-CT (non cancer/testis) subgroup of the melanoma-associated antigen (MAGE) superfamily. The encoded protein is likely associated with apoptosis, cell cycle arrest, growth inhibition or cell differentiation. The


产品名称: MAGH1 Rabbit Polyclonal Antibody
产品货号: APRab13598



protein may be involved in the atRA (all-trans retinoic acid) signaling through the STAT1-alpha (signal transducer and activator of transcription 1-alpha) pathway. [provided by RefSeq, Aug 2013],similarity:Contains 1 MAGE domain.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4°over night

注意事项 (Note)

For research use only .